

## Glucose dehydrogenase (GLD) (NAD(P)-dependent)

Cat no. LDG0023RG

### Product Overview

#### Specification


Appearance	White amorphous powder, lyophilized
Activity	300 U/ mg or more (containing approx. 10% of stabilizers)

#### Properties

Stability	Stable at -20°C for at least one year
Molecular weight	30.3 kDa
Isoelectric point	5.36

#### Assay

##### 1. Assay principle



##### 2. Unit definition

One unit causes the formation of one micromole of NADH per minute under the conditions described below.

##### 3. Reagents

<b>A. Tris-HCl buffer, pH 8.0</b>	0.1 M
<b>B. D-Glucose solution</b>	1.5 M
<b>C. <math>\beta</math>-NAD<sup>+</sup> solution</b>	80 mg/ mL
<b>D. Enzyme diluent</b>	50 mM K-phosphate buffer, pH 7.0 containing 0.1% BSA

#### 4. Procedure

- Prepare the following **working solution** immediately before use and equilibrate at 37°C for approximately 5 minutes (for 4 reactions).

#### Working solution

<b>Tris-HCl buffer, pH 8.0 (Reagent A)</b>	5.2 mL
<b>D-Glucose solution (B)</b>	0.6 mL
<b><math>\beta</math>-NAD<sup>+</sup> solution (C)</b>	0.2 mL
<b>Total</b>	<b>6 mL</b>

- Pipette 1.5 mL of working solution into a tube.
- Add 0.025 mL of the enzyme solution\* and mix by gentle inversion.

Concentration in a reaction	
Tris-HCl buffer	85.25 mM
D-Glucose	147.54 mM
NAD <sup>+</sup>	3.66 mM

- Pipette the mixture into a cuvette (d=1.0 cm).
- Record the increase in optical density at 340 nm against water for 1 to 5 minutes in a spectrophotometer at room temperature and calculate the  $\Delta$ OD per minute from the initial linear portion of the curve ( $\Delta$ OD test). At the same time, measure the blank rate ( $\Delta$ OD blank) by using the same method as the test except that the enzyme diluent is added instead of enzyme solution.

\* Dilute the enzyme in ice-cold enzyme diluent (**Reagent D**) to **0.8–1.2 U/ mL** and store on ice.

(6) Activity can be calculated by using the following formula:

$$\text{Volume activity (U/ mL)} = \frac{\Delta\text{OD/ min } (\Delta\text{OD test} - \Delta\text{OD blank}) \times V_t \times df}{6.22 \times 1.0 \times V_s}$$

$$= \Delta\text{OD/ min} \times 9.807 \times df$$

**Weight activity (U/ mg)** = (U/ mL) × 1/C

Vt: Total volume (1.525 mL)

Vs: Sample volume (0.025 mL)

6.22: Millimolar extinction coefficient of NADH at 340 nm ( $\text{cm}^2$  /micromole)


1.0: Light path length (cm)

df: Dilution factor

C: Enzyme concentration in dissolution (mg/ mL)

### The effect of different conditions on Glucose dehydrogenase


A.



**Figure A. Temperature activity of Glucose dehydrogenase.**

The enzyme reactions in 0.1 M Tris-HCl buffer, pH 8.0, were carried out under different temperature.


B.



**Figure B. pH activity of Glucose dehydrogenase.**

The buffer conditions with various pH values were used in the reaction at 37°C. pH 4.0-6.0, 0.1 M Sodium citrate buffer; pH 7.0, 0.1 M Potassium phosphate buffer; pH 8.0, 0.1 M Tris-HCl buffer; pH 10.0-11.0, 0.1 M Carbonate-bicarbonate buffer.


C.



**Figure C. Thermal stability of Glucose dehydrogenase.**

The enzyme powder was reconstituted by double-distilled water and treated with different temperature for 15 minutes. Final concentration: 23.5 U/ mL

D.



**Figure D. pH stability of Glucose dehydrogenase.** The enzyme powder was reconstituted by double-distilled water and treated with different pH buffer condition for 20 hours. pH 4.0-6.0, 0.1 M Sodium citrate buffer; pH 7.0-8.0, 0.1 M Potassium phosphate buffer; pH 8.5, 0.1 M Tris-HCl buffer; pH 10.0-11.0, 0.1 M Carbonate-bicarbonate buffer.

**Disclaimer**

This product is for research use only and is not intended for diagnostic use.

**LEADGENE BIOMEDICAL, INC.**

No.9, Ln. 147, Zhengbei 1st Rd., Yongkang Dist., Tainan City 710, Taiwan R.O.C. TEL: +886-6-2536677 FAX: +886-6-2531536  
www.leadgenebio.com