

## HyLink™ Aster 594 Labeling Kit, 100 µg\*10 (SpinDesalt Column)

**Catalog Number** LDG0018RC

**Package** 100 µg\*10 / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

The HyLink™ Aster594 Conjugation Kit is tailored for conjugating proteins or antibodies with Aster594, a red fluorescent dye known for its exceptional brightness and photostability. This kit supports small-scale labeling (100 µg) with a straightforward protocol, delivering high efficiency and consistent results. Aster594 is ideal for multicolor fluorescence microscopy, flow cytometry, and advanced imaging applications.

#### Components

Package	Items	Quantity
100 µg × 10	Aster 594	10 vials
	10X Modifier	1 vial
	10X Quencher	1 vial
	SpinDesalt column	10 vials

### Instruction

#### Tainan Headquarter

+886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

---

### Shipping

The product is shipped with polar packs to maintain its stability during transit. Upon receipt, store the SpinDesalt Column at 4°C for optimal performance. All other components should be stored immediately at -20°C or lower for long-term preservation.

---

---

### Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.
- SpinDesalt Column should be store at 2-8°C.

Equilibrate kit to room temperature before use.

---

**Disclaimer :** For Research Use or Further Manufacturing Only.

---

#### Tainan Headquarter

+886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

---

#### Innovation & Research Center

+886-2-27065528

---

#### CLD Center

+886-6-2536677