

HyLink™ Aster 488 Labeling Kit, 100 µg*10 (SpinDesalt Column)

Catalog Number LDG0014RC

Package 100 µg*10 / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

The HyLink™ Aster488 Conjugation Kit enables precise conjugation of proteins or antibodies with Aster488, a bright green fluorescent dye. With optimized protocols for small quantities (100 µg), the kit ensures high-efficiency labeling and strong fluorescence signals suitable for diverse applications, including flow cytometry and imaging. Aster488's superior photostability enhances multicolor analysis, making it a preferred choice for advanced fluorescence-based studies.

Components

Package	Items	Quantity
100 µg × 10	Aster 488	10 vials
	10X Modifier	1 vial
	10X Quencher	1 vial
	SpinDesalt column	10 vials

Instruction

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Shipping

The product is shipped with polar packs to maintain its stability during transit. Upon receipt, store the SpinDesalt Column at 4°C for optimal performance. All other components should be stored immediately at -20°C or lower for long-term preservation.

Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.
- SpinDesalt Column should be store at 2-8°C.

Equilibrate kit to room temperature before use.

Disclaimer : For Research Use or Further Manufacturing Only.