

SARS-CoV-2 Trimeric Spike Protein, His Tag, HEK293

Catalog Number LDG010PHM

Package

100 μg / 1 mg / Customized package

Publications (3)

Specifications

Species of Origin

SARS-CoV-2

Affinity Tag

His Tag (C-term)

Buffer

Liquid. Phosphate buffered saline, pH 7.4.

Molecular weight

The protein has a calculated MW of 137.94 kDa. The protein migrates above 180 kDa under reducing condition (SDS-PAGE analysis).

For full product information, images and publications, please visit our website



please visit our website.

HEK293
Concentration
Please refer to the product label.
Purity
ered saline, pH 7.4.
Possible Section 2007 Participation of 137.94 kDa.
Liquid
Liquid
Liquid

Expression system

Background

Tainan Headquarter

Innovation & Research Center

CLD Center

® +886-6-2536677

& +886-2-<u>27065528</u>

® +886-6-2536677

☑ bd@leadgene.com.tw



Background

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), is a 2019 novel coronavirus, that causes pandemic disease globally. SARS-CoV-2 is a single-stranded RNA virus that belongs to the Coronaviridae family. Structurally, the virus harbors sixteen non-structural proteins (NSP1-NSP16), four structural proteins (spike, membrane, envelope, and nucleocapsid), and nine accessory proteins. The spike protein is a 1273-amino-acids transmembrane. homotrimeric class I fusion protein that contains two functional subunits, receptor-binding subunit S1 (amino acids 14-685) and membrane-fusion subunit S2 (amino acids 686-1273). The trimeric spike protein on the surface of the virion is responsible for initiating viral infection through binding to the angiotensin-converting enzyme 2 (ACE2) receptor on host cells.

Synonyms

Spike glycoprotein, E2, Peplomer protein, S glycoprotein

Sequence Note Ser13-Pro1213

Uniprot ID

#P0DTC2

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -80°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

• -80°C for long-term storage under sterile conditions.

Avoid repeated free-thaw cycles.

Image

Tainan Headquarter

Innovation & Research Center

CLD Center


& +886-6-2536677

® +886-2-27065528

® +886-6-2536677

☑ bd@leadgene.com.tw





SDS-PAGE analysis of recombinant SARS-CoV-2 trimeric spike protein.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

Innovation & Research Center

CLD Center

® +886-6-2536677

& +886-2-2706<u>5528</u>

& +886-6-2536677

 \bowtie bd@leadgene.com.tw