

## HyFect® Transfection Reagent

**Catalog Number** LDG0001RA

**Package** 500 µL / Customized package

“ Publications ( 27 )

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

HyFect® Transfection Reagent is a low cell toxicity and non-liposome reagent optimized for DNA delivery. It provides effectively, reproducibly, and affordable benefits for scientific research.

### Specifications

#### Application

Cell transfection

#### Sample Type

Adherent cells | Suspension cells

### Instruction

#### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

#### Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

**Disclaimer :** For Research Use or Further Manufacturing Only.

**Tainan Headquarter**

+886-6-2536677

bd@leadgene.com.tw

**Innovation & Research Center**

+886-2-27065528

**CLD Center**

+886-6-2536677