

HyFect® Transfection Reagent

Catalog Number LDG0001RA

Package 500 µL / Customized package

“ Publications (27)

For full product information, images and publications, please visit [our website](#).



Overview

Description

HyFect® Transfection Reagent is a low cell toxicity and non-liposome reagent optimized for DNA delivery. It provides effectively, reproducibly, and affordable benefits for scientific research.

Specifications

Application

Cell transfection

Sample Type

Adherent cells | Suspension cells

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677