

## Normal Rabbit IgG

**Catalog Number** LDG0002YD

**Package** 50 µg / 5 mg / 10 mg / Customized package

“ Publications ( 3 )

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Rabbit IgG is isolated and purified from pooled normal rabbit serum by protein A. The purified normal rabbit IgG may be used as an isotype control in various immunoassays, such as ELISA, WB, IP, FACS, ID, and IE.

### Specifications

#### Host

Rabbit

#### Application

ELISA, WB, IP, FACS, ID, IE

#### Concentration

Lot-dependent, please refer to product label.

#### Form

Liquid

#### Clonality

Polyclonal

#### Conjugation

Unconjugated

#### Storage Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

### Instruction

#### Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

---

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

---

### Stability & Storage

"This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C."

---

**Disclaimer :** For Research Use or Further Manufacturing Only.