

Normal Goat IgG

Catalog Number LDG0003YD

Package 5 mg / 10 mg / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Goat IgG is isolated and purified from pooled normal goat serum by protein A. The purified normal goat IgG may be used as an isotype control in various immunoassays, such as ELISA, WB, IP, FACS, ID, and IE.

Specifications

Host

Goat

Clonality

Polyclonal

Application

ELISA, WB, IP, FACS, ID, IE

Conjugation

Unconjugated

Concentration

Lot-dependent, please refer to product label.

Storage Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

Form

Liquid

Instruction

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

☎ +886-6-2536677

✉ bd@leadgene.com.tw

Innovation & Research Center

☎ +886-2-27065528

CLD Center

☎ +886-6-2536677