

Mouse IgG3 Kappa Isotype Control Antibody

Catalog Number LDG0006YC

Package 1mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Mouse IgG3 Kappa isotype control antibody is used to estimate non-specific binding of target primary antibodies.

Specifications

Host

Mouse

Isotype

IgG3

Conjugation

Unconjugated

Storage Buffer

Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

Clonality

Monoclonal

Application

ELISA, WB, IP, FACS, ID, IE

Concentration

1 mg/mL

Form

Liquid

Instruction

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.