


Goat Anti-Rabbit IgG (H+L) Antibody (HRP)

Catalog Number LDG0014YE

Package 1 mL / Customized package

 Publications (32)

For full product information, images and publications, please visit [our website](#).



Overview

Description

Leadgene Goat Anti-Rabbit IgG (H+L) Antibody (HRP) is conjugated with high-activity horseradish peroxidase (HRP).

Product Note

Recommended dilution factor:

ELISA: 1:5000-25000

IHC: 1:500-2500

WB: 1:5000-25000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Goat

Immunogen

Purified rabbit IgG

Application

ELISA, IHC, WB

Clonality

Polyclonal

Reactivity

Rabbit (Minor cross-reactivity to human in higher concentration)

Conjugation

HRP

Storage Buffer

Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300 and 50% glycerol, pH 7.4.

Form

Liquid

Specificity

IgG

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles and protect from light.

Suggestion: Divide antibody into several vials. Keep only vials for usage at 2-8°C.

Disclaimer : For Research Use or Further Manufacturing Only.