

LG-1 2X Diluent Buffer

Catalog Number LDG0017RH

Package 10 mL / 50 mL / 100 mL

“ Publications (1)

For full product information, images and publications, please visit [our website](#).



Overview

Description

LG-1 2X Diluent Buffer is used for dilution of biological samples or antibodies. LG-1 2X Diluent Buffer can reduce matrix effect and non-specific binding in various immunoassays, such as ELISA, WB and IHC.

Specifications

Form

Liquid. LG-1 2X Diluent Buffer was aseptically filtered through a 0.22 μ m filter.

Instruction

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 12 months in lyophilized state from date of receipt.

Disclaimer : For Research Use or Further Manufacturing Only.

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677