

Goat Anti-Rabbit IgG (H+L) Antibody (iFluor 594)

Catalog Number LDG0035YE

Package 500 μ L / 1 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

iFluor 594 dye is a bright, red-fluorescent dye that can be excited by the 561 nm or 594 nm laser. Leadgene Goat Anti-Rabbit IgG (H+L) Antibody (iFluor 594) is conjugated with high-activity iFluor 594.

Product Note

Recommended dilution factor:

IFA: 1:200-1000

FACS: 1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Goat

Immunogen

Purified rabbit IgG

Application

IFA, FACS

Clonality

Polyclonal

Reactivity

Rabbit (Minor cross-reactivity to human in higher concentration)

Conjugation

iFluor 594

Concentration

0.5 mg/mL

Storage BufferPhosphate Buffered Saline containing 1% BSA,
0.03% ProClin 300, pH 7.4.**Specificity**

IgG

Form

Liquid

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 6 months under sterile conditions from date of receipt.

Do not freeze and protect from light.

Disclaimer : For Research Use or Further Manufacturing Only.