

## Goat Anti-Rat IgG (H+L) Antibody (FAM)

<b>Catalog Number</b>	LDG0045YE
<b>Package</b>	500 µL / 1 mL / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

FAM is a fluorescent dye that can be excited by wavelength of 495 nm. Leadgene Goat Anti-Rat IgG (H+L) Antibody (FAM) is conjugated with high-activity FAM.

#### Product Note

Recommended dilution factor:  
IFA: 1:200-1000  
FACS: 1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

### Specifications

#### Host

Goat

#### Immunogen

Purified rat IgG

#### Application

IFA, FACS

#### Clonality

Polyclonal

#### Reactivity

Rat

#### Conjugation

5(6)-FAM

**Concentration**

0.5 mg/mL

**Specificity**

IgG

**Storage Buffer**

Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.

**Form**

Liquid

**Instruction****Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

**Stability & Storage**

This product is stable after storage at:

- 2-8°C for 6 months under sterile conditions from date of receipt.

Do not freeze and protect from light.

**Disclaimer :** For Research Use or Further Manufacturing Only.