

Donkey Anti-Goat IgG (H+L) Antibody (TAMRA)

Catalog Number LDG0048YE

Package 500 µL / 1 mL / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

TAMARA is a fluorescent dye for DNA, peptides and proteins. Leadgene Donkey Anti-Goat IgG (H+L) Antibody (TAMRA) is conjugated with high-activity TAMARA

Product Note

Recommended dilution factor:

IFA: 1:200-1000

FACS: 1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Donkey

Immunogen

Purified goat IgG

Application

IFA, FACS

Clonality

Polyclonal

Reactivity

Goat (Minor cross-reactivity to human in higher concentration)

Conjugation

5(6)-TAMRA

Concentration

0.5 mg/mL

Storage BufferPhosphate Buffered Saline containing 1% BSA,
0.03% ProClin 300, pH 7.4.**Specificity**

IgG

Form

Liquid

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 6 months under sterile conditions from date of receipt.

Do not freeze and protect from light.

Disclaimer : For Research Use or Further Manufacturing Only.