

Goat Anti-Chicken IgY Antibody

Catalog Number LDG0052YE

Package 1 mg / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Goat anti-chicken IgY antibody is purified from goats and specifically reacts with chicken IgY antibodies.

Specifications

Host

Goat

Immunogen

Purified chicken IgY

Concentration

Lot-dependent, please refer to product label.

Specificity

IgY

Clonality

Polyclonal

Conjugation

Unconjugated

Storage Buffer

Phosphate Buffered Saline containing 0.03% ProClin300, pH 7.4.

Form

Liquid

Instruction

Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

Innovation & Research Center

+886-2-27065528

CLD Center

+886-6-2536677

Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage

This product is stable after storage at:


- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.


Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.

Image



ELISA titration of goat anti-chicken IgY antibody



SDS-PAGE of purified antibody under reduced condition.

Lane M: Protein marker

Lane 1: Goat anti-Chicken IgY

(H+L) Antibody (2 µg, reduced)

Disclaimer : For Research Use or Further Manufacturing Only.