

Anti-MPXV A29L Antibody [Clone 19-1]

Catalog Number LDG0189YA

Package 100 µg / Customized package

For full product information, images and publications, please visit [our website](#).



Overview

Description

Monkeypox is a zoonotic disease caused by the Monkeypox virus (MPXV), a DNA virus that belongs to the Orthopoxvirus genus, which also includes Variola virus (smallpox) and Cowpox virus. The MPXV A29L protein is a structural protein and is homologous to the vaccinia virus Copenhagen A27 protein. It is a critical component of the virus particle and is essential for viral entry by interacting with heparin on the cell surfaces. Anti-MPXV A29L Antibody [Clone 19-1] recognize Monkeypox virus (MPXV) A29L protein.

Product Note

Recommended dilution factor:

ELISA: 1:5000-20000

IFA: 1:500-1:1000

FACS: Assay dependent

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Mouse

Clonality

Monoclonal

Isotype

IgG1

Clone Name

clone 19-1

Immunogen

Recombinant MPXV A29L protein

Reactivity

Monkeypox virus

Application

ELISA, WB, IFA

Concentration

1 mg/mL

Specificity

A29L protein

Conjugation

Unconjugated

Storage Buffer

Phosphate Buffered Saline containing 0.03% ProClin 300, pH 7.4.

Form

Liquid

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

Stability & Storage




This product is stable after storage at:


- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.


Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.

Image**Tainan Headquarter** +886-6-2536677 bd@leadgene.com.tw**Innovation & Research Center** +886-2-27065528**CLD Center** +886-6-2536677



Immunofluorescence analysis of Anti-MPXV A29L Antibody [Clone 19-1].



Western blotting analysis of Anti-MPXV A29L Antibody [Clone 19-1]
Western blot: A29L expressed 293T cell lysates were stained with mouse anti-MPXV A29L monoclonal antibody. Lane 1: 293T mock cell lysate, Lane 2: 293T MPXV cell lysate.

Disclaimer : For Research Use or Further Manufacturing Only.