

## Goat Anti-Human IgG (Fc Fragment) Antibody

**Catalog Number** LDG0011YE

**Package** 1 mg / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Leadgene Goat anti-Human IgG, Fc Fragment is purified from goat antiserum by immunoaffinity chromatography. Minimal cross-reactivity with other isotype or species.

### Specifications

#### Host

Goat

#### Clonality

Polyclonal

#### Immunogen

Gamma immunoglobulin Fc region

#### Reactivity

Human

#### Conjugation

Unconjugated

#### Concentration

1 mg/mL

#### Buffer

Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.

#### Specificity

IgG, Fc Fragment

#### Form

Liquid

### Instruction

#### Tainan Headquarter

+886-6-2536677

bd@leadgene.com.tw

#### Innovation & Research Center

+886-2-27065528

#### CLD Center

+886-6-2536677

---

### Shipping

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

---

### Stability & Storage

This product is stable after storage at:

- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.

Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.

---

**Disclaimer :** For Research Use or Further Manufacturing Only.

---

#### Tainan Headquarter

☎ +886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

---

#### Innovation & Research Center

☎ +886-2-27065528

---

#### CLD Center

☎ +886-6-2536677