

## Anti-mCherry Antibody Sepharose

**Catalog Number** LDG0021RD

**Package** 0.25 mL / 2 mL / 5 mL / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Anti-mCherry VHH antibody was conjugated to NHS-sepharose beads. Leadgene Mouse anti-mcherry mAb sepharose has 6-10 mg mAb conjugated per mL of settled sepharose.

#### Product Note

The optimal protein concentration of the cell lysate should be determined by the end user.

### Specifications

#### Clonality

VHH

#### Isotype

IgG2a

#### Immunogen

Recombinant mCherry protein

#### Application

Immunoprecipitation, Immunoaffinity Purification

#### Conjugation

Sepharose

#### Specificity

Recognizes RFP, DsRed and mCherry tag on either the N- or C-terminus of the recombinant fusion protein.

#### Sample Type

Cell lysate

#### Form

Gel Immobilized 50% slurry in 20% ethanol.

---

**Binding capacity**

The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of RFP, mCherry or DsRed fusion protein.

---

**Instruction**

---

**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C.

---

---

**Stability & Storage**

This product is stable after storage at:

- 2-8°C for 12 months under sterile conditions from date of receipt.

Do not freeze.

---

**Disclaimer :** For Research Use or Further Manufacturing Only.

---

**Tainan Headquarter**

☎ +886-6-2536677

✉ [bd@leadgene.com.tw](mailto:bd@leadgene.com.tw)

---

**Innovation & Research Center**

☎ +886-2-27065528

---

**CLD Center**

☎ +886-6-2536677