

## Goat Anti-Mouse IgG (H+L) Antibody (iFluor 594)

**Catalog Number** LDG0034YE

**Package** 500  $\mu$ L / 1 mL / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

iFluor 594 dye is a bright, red-fluorescent dye that can be excited by the 561 nm or 594 nm laser. Leadgene Goat Anti-Mouse IgG (H+L) Antibody (iFluor 594) is conjugated with high-activity iFluor 594.

#### Product Note

Recommended dilution factor:

IFA: 1:200-1000

FACS: 1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

### Specifications

#### Host

Goat

#### Immunogen

Purified mouse IgG

#### Application

IFA, FACS

#### Clonality

Polyclonal

#### Reactivity

Mouse (Minor cross-reactivity to rat)

#### Conjugation

iFluor 594

**Concentration**

0.5 mg/mL

**Buffer**Phosphate Buffered Saline containing 1% BSA,  
0.03% ProClin 300, pH 7.4.**Specificity**

IgG

**Form**

Liquid

**Instruction****Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

**Stability & Storage**

This product is stable after storage at:

- 2-8°C for 6 months under sterile conditions from date of receipt.

Do not freeze and protect from light.

**Disclaimer :** For Research Use or Further Manufacturing Only.