

Goat Anti-Mouse IgG (H+L) Antibody (FAM)

Catalog Number LDG0040YE

Package 500 μ L / 1 mL / Customized package

“ Publications (1)

For full product information, images and publications, please visit [our website](#).



Overview

Description

FAM is a fluorescent dye that can be excited by wavelength of 495 nm. Leadgene Goat Anti-Mouse IgG (H+L) Antibody (FAM) is conjugated with high-activity FAM.

Product Note

Recommended dilution factor:

IFA: 1:200-1000

FACS: 1:200-1000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

Specifications

Host

Goat

Immunogen

Purified mouse IgG

Application

IFA, FACS

Clonality

Polyclonal

Reactivity

Mouse (Minor cross-reactivity to rat)

Conjugation

5(6)-FAM

Concentration

0.5 mg/mL

Specificity

IgG

Buffer

Phosphate Buffered Saline containing 1% BSA, 0.03% ProClin 300, pH 7.4.

Form

Liquid

Instruction**Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at 2-8°C for long term storage.

Stability & Storage

This product is stable after storage at:

- 2-8°C for 6 months under sterile conditions from date of receipt.

Do not freeze and protect from light.

Disclaimer : For Research Use or Further Manufacturing Only.