

## Anti-Influenza A Virus NP Antibody [Clone IA02]

**Catalog Number** LDG0144YA

**Package** 100 µg / Customized package

For full product information, images and publications, please visit [our website](#).



### Overview

#### Description

Anti-Influenza A virus NP Antibody only recognizes nucleocapsid protein (NP) of influenza A viruses but not influenza B viruses. Influenza A viruses are RNA viruses, and their subtypes are labeled according to an H number (for the type of hemagglutinin) and an N number (for the type of neuraminidase).

#### Product Note

Recognize Influenza A virus NP in Lateral Flow, when clone IA02 antibody was paired with Anti-Influenza A virus NP Antibody [clone IA01] (cat. LDG0143YA).

Recommended dilution factor:  
ELISA: 1:5000-20000

Note: Working dilution for specific application should be determined by the investigator to obtain the best conditions.

### Specifications

#### Host

Mouse

#### Clonality

Monoclonal

#### Isotype

IgG1, IgG2a

#### Clone Name

clone IA02

**Immunogen**

Recombinant Influenza A virus nucleocapsid protein

**Application**

ELISA, IFA, CLIA, LFIA

**Concentration**

1 mg/mL

**Specificity**

Nucleocapsid protein

**Reactivity**

Influenza A virus

**Conjugation**

Unconjugated

**Buffer**

Phosphate Buffered Saline containing 0.2% ProClin 300, pH 7.4.

**Form**

Liquid

**Instruction****Shipping**

The product is shipped with polar packs. Upon receipt, store it immediately at -20°C or lower for long term storage.

**Stability & Storage**

This product is stable after storage at:


- 2-8°C for 2 weeks under sterile conditions from date of receipt.
- -20°C or -80°C for 12 months under sterile conditions from date of receipt.

Avoid repeated freeze/thaw cycles.


Suggestion: Divide antibody into several vials.

Keep only vials for usage at 2-8°C.

**Image**



Immunofluorescence analysis of Anti-Influenza A virus NP Antibody [clone IA02] (1:500)



ELISA assay of Anti-Influenza A virus NP Antibody [clone IA02]

**Disclaimer :** For Research Use or Further Manufacturing Only.